

## Surge 15W – 30W UV-C Lamb

Air and Surface Disinfection Device / Hava ve Yüzey Dezenfeksiyon Cihazı

### Datasheet / Teknik Özellikler



### FIELDS of USAGE

It is designed only for **hospitals, clinics** and **pharmacies**. It should not be used anywhere else. Warnings should be observed while using. **UV-C light** sources can **cause permanent damage to human health**.

MODEL	Surge 15W – 30W UV-C Lambs
TYPE	Wall Type Surface and Air Disinfection Device
CLASS	-
PHYSICAL AND TECHNICAL FEATURES	
EXT. DIMENSIONS (W-L-H) (mm)	121x535x126 (15 W) 121x990x126 (30 W)
DEVICE WIGHT (kg)	2.7 (15 W) 4.1 (30 W)
LAMB BASE TYPE	G13
LAMB TYPE	TL-D
BALAST TYPE	Electronic
LIGHT FLUX (lm)	960 (15 W) 2175 (30 W)

INPUT CURRENT (A <sub>RMS</sub> ) @230 V <sub>RMS</sub>	0.07 (15 W) 0.14 (30 W)
UV-C LAMB FEATURES	1 Adet 15W <b>UV-C</b> 253.7 nm (15 W) 1 Adet 30W <b>UV-C</b> 253.7 nm (30 W)
ELECTRICAL FEATURES	220-240V <sub>RMS</sub> 50 Hz
PROTECTION CLASS	IP20
<b>CERTIFICATIONS</b>	
CE CERTIFICATE	Yes
EN-60598-2-1	Yes
<b>INSTALLATION FEATURE</b>	
Thanks to its hanging apparatus, it can be attached to the wall.	
<b>WARNINGS</b>	
Ultraviolet lamps <b><u>should not be directly looked at or touched</u></b> . The front shield of UV-C lamps should be kept closed. It should be mounted on 120 cm and above. It should <b><u>never be operated</u></b> in environments with <b>people or living things</b> .	
<b>CARE PERIOD AND CARE KIT</b>	
UV LAMB	It should be cleaned weekly with lint-free microfiber cloth and alcohol while the light is off. It should be changed once in 1 year.